

2019-nCoV Ag Saliva Rapid Test Card (Immunochromatography)

Catalog Number:

0589C4X001 0589C4X005 0589C4X010 0589C4X015 0589C4X020

INTENDED USE

The Test Card is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigen from 2019-nCoV in saliva specimens directly collected from individuals who are suspected of COVID-19 by their healthcare provider within the first 7 days of symptom onset.

Results are for the identification of 2019-nCoV nucleocapsid protein antigen. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results should be treated as presumptive, and do not rule out 2019-nCoV infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management.

SUMMARY AND EXPLANATION

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

PRINCIPLE OF THE TEST

This Card uses double-antibody sandwich to legally detect the antigen of novel coronavirus (2019-nCoV) in saliva samples. During detection, the gold labeled anti-2019-nCoV monoclonal antibody in the labeling pad binds to the 2019-nCoV antigen in the sample to form a complex, and the reaction complex moves forward along the nitrocellulose membrane under the action of chromatography, it is captured by the anti-2019-nCoV monoclonal antibody pre-coated by the detection zone (T) on the nitrocellulose membrane, and finally a red color reaction line is formed in the T zone. If the sample does not contain 2019-nCoV antigen, a red color reaction line cannot be formed in the T zone. Regardless of whether the sample to be tested contains 2019-nCoV antigen, a red reaction line will always form in the quality control area (C).

MATERIALS AND COMPONENTS

Materials provided with the test kits

| Specifications | 0589C 4X001 | 0589C 4X005 | 0589C 4X010 | 0589C 4X015 | 0589C 4X020 |
|------------------------------|----------------|----------------|----------------|----------------|----------------|
| Ingredients | | | | | |
| Test Card | 1 | 5 | 10 | 15 | 20 |
| Saliva Swab | 1 | 5 | 10 | 15 | 20 |
| Instructions for use | 1 | 1 | 1 | 1 | 1 |
| Quick Reference Instructions | NA | 1 | 1 | 1 | 1 |

Materials required but not provided

1. Timer

STORAGE AND STABILITY

1. Store the test card as packaged between 2-30°C.
2. The Test Card is stable until the expiration date printed on the outer packing, the product will be expired after 24 months.
3. Do not use beyond the expiration date.
4. Do not freeze any contents of the test.
5. The test card must remain in the sealed pouch until use.

SAMPLE REQUIREMENTS


1. Insert the sponge end of the saliva swab into mouth. Actively swab the inside of the mouth and tongue to collect oral fluid.
2. Remove the saliva swab from the mouth when the sponge fill with saliva and become soft, or the indicator turns blue.
3. Do not eat, drink or smoke prior to the test for at least 30 Minutes.
4. The samples should be used as soon as possible after collected.
5. Samples should not be inactivated.

NOTE:

*When sampling, gently hold it in mouth and let saliva naturally adsorb on the sponge.

*Don't bite the sponge with teeth.


*Any saliva specimen is appropriate for testing but the saliva specimen collected in the morning, before mouth rinsed, eating or drinking, is recommended.



TEST PROCEDURE

Before test, please read the instructions carefully.

1. Take the test card to equilibrate to room temperature.
2. Unpack the aluminum foil bag, place the test card horizontally on the table and mark it.
3. Insert the saliva swab into the test card holder and push down saliva swab. The bump at the end of the saliva swab must be into the hole of the test card holder.
4. As the test begins to work, the purple color move across the result window in the center of the test device.
5. Wait for 10 minutes and read the results.



INTERPRETATION OF TEST RESULTS

This product can only perform qualitative analysis on the detection object.

Positive Result:


If both C and T lines are visible within 10 minutes, the test result is positive and valid.

Negative Result:

If test area (T line) has no color and the control area displays a colored line, the result is negative and valid.

Invalid Result:

The test result is invalid if a colored line does not form in the control region. The sample must be re-tested, using a new test card.



LIMITATIONS

1. The result of the test card should not be taken as a confirmed diagnosis, for clinical reference only. Judgement should be made along with RT-PCR results, clinical symptoms, epidemiological information and further clinical data.
2. Test Card performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
3. The test card must be equilibrated to room temperature (18°C~26°C) before used, otherwise the results may be incorrect
4. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.
5. Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
6. React less than 10 minutes may lead a false negative result; React more than 10 minutes may lead a false positive result.
7. Positive test results do not rule out co-infections with other pathogens.
8. Negative test results are not intended to rule in other viral or bacterial infections.
9. Negative results should be treated as presumptive and confirmed with a molecular assay.
10. Clinical performance was evaluated with fresh samples.
11. Users should test specimens as quickly as possible after specimen collection.

PERFORMANCE CHARACTERISTIC

1. Clinical Verification

The performance of Test Card was established with 243 sample collected from symptomatic patients, who with symptoms onset within 7 days.

| 2019-nCoV Ag Saliva Rapid Test Card (Immunochromatography) | Comparative RT-PCR Test Result | | |
|--|--------------------------------|--------------|-------|
| | Positive (+) | Negative (-) | Total |
| Detected Positive | 110 | 2 | 112 |
| Detected Negative | 5 | 126 | 131 |
| Total | 115 | 128 | 243 |
| Sensitivity | 95.65%, 95% CI (90.22,98.13) | | |
| Specificity | 98.44%, 95% CI (94.48, 99.57). | | |
| Accuracy | 97.12%, 95% CI (94.17,98.60) | | |

The performance of Test Card with positive results stratified by the comparative method cycle threshold (Ct) counts were collected and assessed to better understand the correlation of assay performance to the cycle threshold. As presented in the table below, the positive agreement of the Test Card is higher with samples of a Ct count <25.

| 2019-nCoV Ag Saliva Rapid Test Card (Immunochromatography) | Comparative RT-PCR Method (Positive by Ct Value) | |
|--|--|------------------|
| | Positive>25 (Ct<=25) | Positive (Ct>25) |
| Detected Positive | 91 | 19 |
| Total | 92 | 23 |
| Positive agreement | 98.91% | 82.60% |

2. Limit of Detection

The experimental results show that for the virus culture concentration above 100 TCID₅₀/mL, the positive rate of detection is greater than or equal to 95%. For the virus culture concentration of 50 TCID₅₀/mL and below, the positive rate of detection is lower than 95%. So, the limit of detection of the Test Card is 100 TCID₅₀/mL.

3. Cross-reactivity

Cross-reactivity of the test card was evaluated. The results showed no cross reactivity with the following specimen.

| No. | Specimen type | Conc. |
|-----|-------------------------------|--|
| 1 | HCoV-HKU1 | 10 ⁵ TCID ₅₀ /mL |
| 2 | Staphylococcus aureus | 10 ⁶ CFU / mL |
| 3 | Streptococcus pyogenes | 10 ⁶ CFU / mL |
| 4 | Measles virus | 10 ⁵ TCID ₅₀ /mL |
| 5 | Paramyxovirus parotitis | 10 ⁵ TCID ₅₀ /mL |
| 6 | Adenovirus 3 | 10 ⁵ TCID ₅₀ /mL |
| 7 | Mycoplasma pneumoniae | 10 ⁶ CFU / mL |
| 8 | Parainfluenza virus 2 | 10 ⁵ TCID ₅₀ /mL |
| 9 | Human Metapneumovirus (hMPV) | 10 ⁵ TCID ₅₀ /mL |
| 10 | Human coronavirus OC43 | 10 ⁵ TCID ₅₀ /mL |
| 11 | Human coronavirus 229E | 10 ⁵ TCID ₅₀ /mL |
| 12 | Human coronavirus NL63 | 10 ⁴ TCID ₅₀ /mL |
| 13 | MERS-Coronavirus EMC/2012 | 10 ⁴ TCID ₅₀ /mL |
| 14 | Bordetella parapertussia | 10 ⁶ CFU / mL |
| 15 | Influenza B (Victoria strain) | 10 ⁵ TCID ₅₀ /mL |
| 16 | Influenza B (Y strain) | 10 ⁵ TCID ₅₀ /mL |
| 17 | Influenza A (H1N1 2009) | 10 ⁵ TCID ₅₀ /mL |
| 18 | Influenza A (H3N2) | 10 ⁵ TCID ₅₀ /mL |
| 19 | Avian influenza virus (H7N9) | 10 ⁵ TCID ₅₀ /mL |
| 20 | Avian influenza virus (H5N1) | 10 ⁵ TCID ₅₀ /mL |
| 21 | Epstein-Barr virus | 10 ⁵ TCID ₅₀ /mL |

| | | |
|----|-----------------------------|--|
| 22 | Enterovirus CA16 | 10 ⁵ TCID ₅₀ /mL |
| 23 | Rhinovirus | 10 ⁵ TCID ₅₀ /mL |
| 24 | Respiratory syncytial virus | 10 ⁵ TCID ₅₀ /mL |
| 25 | Streptococcus pneumoniae | 10 ⁶ CFU / mL |
| 26 | Candida albicans | 10 ⁶ CFU / mL |
| 27 | Chlamydia pneumoniae | 10 ⁶ CFU / mL |
| 28 | Bordetella pertussis | 10 ⁶ CFU / mL |
| 29 | Pneumocystis jirovecii | 10 ⁶ CFU / mL |
| 30 | Mycobacterium tuberculosis | 10 ⁶ CFU / mL |
| 31 | Legionella pneumophila | 10 ⁶ CFU / mL |

4. Interference Substances

The test results do not be interfered with the substance at the following concentration:

| No. | Interference substances | Conc. |
|-----|--------------------------|----------|
| 1 | Whole Blood | 4% |
| 2 | Ibuprofen | 1mg / mL |
| 3 | Tetracycline | 3µg / mL |
| 4 | Chloramphenicol | 3µg / mL |
| 5 | Erythromycin | 3µg / mL |
| 6 | Tobramycin | 5% |
| 7 | Throat spray (Menthol) | 15% |
| 8 | Mupirocin | 10mg/mL |
| 9 | Throat lozenge (Menthol) | 1.5mg/mL |
| 10 | Tamiflu (Oseltamivir) | 5mg/mL |

| | | |
|----|---|----------|
| 11 | Naphthoxoline hydrochloride nasal drops | 15% |
| 12 | Mucin | 0.50% |
| 13 | Fisherman's Friend | 1.5mg/mL |
| 14 | Compound Benzocain Gel | 1.5mg/mL |
| 15 | Cromoglycate | 15% |
| 16 | Sinex (Phenylephrine Hydrochloride) | 15% |
| 17 | Afrin (Oxymetazoline) | 15% |
| 18 | Fluticasone propionate spray | 15% |

5. Precision

- Test 10 replicates of negative and positive by using the reference materials of enterprises. The negative agreement and the positive agreement were 100%.
- Test three different lots kits including positive and negative reference materials of enterprises. The negative results and the positive results were 100%.

















6. Hook Effect

The Test Card was tested up to 1.6 × 10⁶ TCID₅₀/ml of heat-inactivated 2019-nCoV strain and no high-dose effect was observed.


PRECAUTIONS

- For in vitro diagnostic use.
- Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used test contents.
- Use of Nitrile, Latex (or equivalent) gloves is recommended when handling patient samples.
- Do not reuse the used Test Card or saliva swab.
- Should never open the foil pouch of the Test Card exposing it to the ambient environment until the Test Card is ready for immediate use.
- Discard and do not use any damaged or dropped Test Card or material.
- Inadequate or inappropriate sample collection, storage, and transport may yield false test results.
- Sample collection and handling procedures require specific training and guidance.
- To obtain accurate results, do not use visually bloody or overly viscous samples.
- To obtain accurate results, an opened and exposed Test Card should not be used.
- Testing should be performed in an area with adequate ventilation.
- Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this test.
- Wash hands thoroughly after handling.

KEY TO SYMBOLS USED

| | | | |
|---|------------------------------|--|------------------------------------|
|  COMPONENT | Materials Included |  SWAB | Saliva Swab |
|  CASSETTE | Test Card |  | Date of Manufacturer |
|  IFU | Instructions for Use |  | Do Not Reuse |
|  | Consult Instructions For Use |  REF | Catalogue Number |
|  | Store at 2°C~30°C |  | Keep away from Sunlight |
|  | Expiration Date |  | Tests per Kit |
|  | Manufacturer |  IVD | In Vitro Diagnostic Medical Device |
|  LOT | Lot Number |  | Keep Dry |

 **Guangzhou Decheng Biotechnology Co., LTD**
Room 218, Building 2, No.68, Nanxiang Road, Science City, Huangpu District, 510000, Guangzhou P.R.China
TEL:+86-020-82557192
service@dochebio.com
www.dochebio.com

 Caretechion GmbH
Niederrheinstr. 71, 40474
Duesseldorf, Germany.



Doc No.: DC-IN-0589C01 Ver 1.6 EN
Rel.:2021/02/20